


Sequencing Technologies

BIOL 432


Moore's Law

Microprocessor Transistor Counts 1971-2011 & Moore's Law




Moore's Law vs. Sequencing Technology

Cost per Raw Megabase of DNA Sequence



Moore's Law vs. Sequencing Technology

Cost per Raw Megabase of DNA Sequence



Sequencing Overview – Match platform with description

Sanger

First-generation sequencing

Illumina

Second-generation sequencing

PacBio

Third-generation sequencing

Nanopore

Next-generation sequencing

(Also: IonTorrent, 454, SOLID)

Sequencing-by-synthesis

Nanopore sequencing

High-throughput sequencing

Short-read sequencing


Long-read sequencing

Sanger Method

First-generation sequencing


Review: Polymerase chain reaction

Polymerase chain reaction - PCR




1st Generation Sequencing: Sanger Method

= 'Sequencing-by-synthesis'



Beckman Coulter CEQ 8000

Sequencing Read Archive (NCBI)





Blue and green =
'Next Generation Sequencing'

Second-generation sequencing


Short-read, high-throughput sequencing

Next-Generation Sequencing: Typical Workflow



Sequencing Library Preparation




1. Extract & purify DNA
2. Fragment to target size (75-750 bp)
3. Strand isolation
4. Clonal Amplification
5. Nucleotide detection





Fragment Sizes (partial list)

Platform	Instrument	Mreads	Length (bp)	Gbp	Type
Illumina	NovaSeq 6000 S4	10000	300	3000	SR & PE
Illumina	NextSeq 500 High-Output	400	300	120	SR & PE
Illumina	HiSeq X	375	300	112.5	PE
Illumina	HiSeq High-Output v4	250	250	62.5	SR & PE
Illumina	MiSeq v3	25	600	15	SR & PE
Illumina	MiniSeq High-Output	25	300	7.5	SR & PE
Ion	Proton I	60	200	12	SR
Ion	PGM 318	4	400	1.6	SR
Ion	PGM 316	2	400	0.8	SR
Ion	PGM 314	0.4	400	0.16	SR
PacBio	PacBio Sequel	0.37	20000	7.4	SR
PacBio	PacBio RS II (P6)	0.055	15000	0.825	SR
Roche 454	GS FLX+ / FLX	0.7	700	0.49	SR
SOLiD	5500xl W	267	100	26.7	SR & PE



Sequence isolation (and cloning)




Ion Torrent



2nd Gen: 454 Sequencing (Roche; deprecated)




2nd Gen: 454 Sequencing (Roche; deprecated)



CSB2008 August 2008


Rothberg & Leomon 2008

SOLiD Sequencing (ABI)




a


SOLiD™ substrate



DI base probes




1. Prime and ligate




2. Image



3. Cap unextended strands



4. Cleave off fluor




8. Repeat Reset with , n-2, n-3, n-4 primers


5. Repeat steps 1-4 to extend sequence





6. Primer reset




7. Repeat steps 1-5 with new primer



Illumina flow cells (a micro-array technology)




Illumina sequencing (formerly Solexa)






1. Prepare genomic DNA
2. Attach DNA to surface
3. Bridge amplification
4. Fragment become double stranded
5. Denature the double stranded molecules
6. Complete amplification

Illumina sequencing



7. Determine first base
8. Image first base
9. Determine second base
10. Image second base
11. Sequence reads over multiple cycles
12. Align data

Illumina Devices (benchtop)

	 iSeq 100	 MiniSeq	 MiSeq Series	 NextSeq 550 Series	 NextSeq 1000 & 2000
Popular Applications & Methods	Key Application	Key Application	Key Application	Key Application	Key Application
Large Whole-Genome Sequencing (human, plant, animal)					
Small Whole-Genome Sequencing (microbe, virus)	●	●	●	●	●
Exome & Large Panel Sequencing (enrichment-based)				●	●
Targeted Gene Sequencing (amplicon-based, gene panel)	●	●	●	●	●
Single-Cell Profiling (scRNA-Seq, scDNA-Seq, oligo tagging assays)				●	●
Transcriptome Sequencing (total RNA-Seq, mRNA-Seq, gene expression profiling)				●	●
Targeted Gene Expression Profiling	●	●	●	●	●
miRNA & Small RNA Analysis	●	●	●	●	●
DNA-Protein Interaction Analysis (ChIP-Seq)			●	●	●
Methylation Sequencing				●	●
16S Metagenomic Sequencing		●	●	●	●
Metagenomic Profiling (shotgun metagenomics, metatranscriptomics)				●	●
Cell-Free Sequencing & Liquid Biopsy Analysis				●	●

Illumina Devices (industrial scale)



NextSeq 550 Series +



NextSeq 1000 & 2000



NovaSeq 6000

Popular Applications & Methods	Key Application	Key Application	Key Application
Large Whole-Genome Sequencing (human, plant, animal)			
Small Whole-Genome Sequencing (microbe, virus)			
Exome & Large Panel Sequencing (enrichment-based)			
Targeted Gene Sequencing (amplicon-based, gene panel)			
Single-Cell Profiling (scRNA-Seq, scDNA-Seq, oligo tagging assays)			
Transcriptome Sequencing (total RNA-Seq, mRNA-Seq, gene expression profiling)			
Chromatin Analysis (ATAC-Seq, ChIP-Seq)			
Methylation Sequencing			
Metagenomic Profiling (shotgun metagenomics, metatranscriptomics)			
Cell-Free Sequencing & Liquid Biopsy Analysis			

Third-generation sequencing

Long-read sequencing

3rd Generation sequencing: SMRT

SMRT = single-molecule real-time sequencing

Start with high-quality double stranded DNA



Ligate SMRTbell adapters and size select




Anneal primers and bind DNA polymerase





Circularized DNA is sequenced in repeated passes

The polymerase reads are trimmed of adapters to yield subreads


Consensus is called from subreads



	Sequel IIe System	Sequel II System	Sequel System
Supported SMRT Cell	SMRT Cell 8M	SMRT Cell 8M	SMRT Cell 1M
Number of HiFi Reads >99%* Accuracy	Up to 4,000,000	Up to 4,000,000	Up to 500,000
Sequencing Run Time per SMRT Cell	Up to 30 hrs	Up to 30 hrs	Up to 20 hrs



3rd Generation sequencing: Nanopore



Nanopore Sequencing Comparison (2019)

Platform	Instrument	Mreads	Length (bp)	Gbp	Type
Illumina	NovaSeq 6000 S4	10000	300	3000	SR & PE
Illumina	NextSeq 500 High-Output	400	300	120	SR & PE
Illumina	HiSeq X	375	300	112.5	PE
Illumina	HiSeq High-Output v4	250	250	62.5	SR & PE
Illumina	MiSeq v3	25	600	15	SR & PE
Illumina	MiniSeq High-Output	25	300	7.5	SR & PE
Oxford Nanopore	MinION	1M+	20 SR		
Oxford Nanopore	PromethION			1000	SR
Ion	Proton I	60	200	12	SR
Ion	PGM 318	4	400	1.6	SR
Ion	PGM 316	2	400	0.8	SR
Ion	PGM 314	0.4	400	0.16	SR
PacBio	PacBio Sequel	0.37	20000	7.4	SR
PacBio	PacBio RS II (P6)	0.055	15000	0.825	SR
Roche 454	GS FLX+ / FLX	0.7	700	0.49	SR
SOLiD	5500xl W	267	100	26.7	SR & PE

Review

```
A set of bioinformatics packages for R  
install.packages("BiocManager")  
Update all/some/none [a/s/n] choose n  
library(BiocManager)  
install(c("sangerseqR","annotate"))
```

Working in groups (15 mins):

Stretch and divide into working groups

Summarize sequencing-by-synthesis (SBS) with Illumina

Review key concepts:

1. How a flow cell works
2. Contrast Sanger with SBS sequencing

Try flowcharts or cartoons to simplify & summarize

BRAINSTORM:

What are the main benefits & limitations of each technology?

Why is coding valuable for 2nd generation sequencing?